

# Cleaning and Disinfecting the ComfortGel™ Mask




## CLEANING

Clean the mask prior to using any of the disinfection or sterilization methods listed below.


- Remove the headgear and spacer, if appropriate.
- Soak the parts in a commercially available enzymatic cleaner (e.g., Klenzyme™).
- Clean the mask with a soft bristle brush in a solution of cool tap water and a commercially available anionic detergent (e.g., Manu Klenz™). Do not use cleaning products that contain conditioners or moisturizers because they will leave a residue.
- Rinse thoroughly under cool running tap water. Air dry.
- If adhering substances have not been adequately removed, replace the mask.

## DISINFECTION/STERILIZATION PROCESS

### PASTEURIZATION

RESPIRONICS MASK	TEMPERATURE	EXPOSURE TIME	RESULT	NO. CYCLES TESTED	NOTES
ComfortGel 	70° C	30 minutes	High level disinfection	30	<ul style="list-style-type: none"> <li>• Inspect the mask after processing. If any components are damaged, replace the mask.</li> <li>• Tested for durability at 75° C.</li> </ul>

### Cidex

RESPIRONICS MASK	EXPOSURE TIME	RINSE	RESULT	NO. CYCLES TESTED	NOTES
ComfortGel 	45 minutes	Per manufacturer's instruction	High level disinfection	30	<ul style="list-style-type: none"> <li>• Inspect the mask after processing. If any components are damaged, replace the mask.</li> <li>• Slight discoloration of the mask cushion after processing is normal.</li> </ul>

Contact Customer Service to obtain mask disinfecting instructions for other Respironics masks (#1011138).