

REMstar[®] Auto Smart CPAP System (P/N: 1007381)

SPECIFICATIONS

Pressure Range	4 to 20 cm H ₂ O
Ramp Time (CPAP Mode Only)	0 to 45 min. (5-min. increments)
Starting Ramp Pressure (CPAP Mode Only)	4 to CPAP pressure; patient-adjustable
Dimensions	9.45" (L) x 6.69" (W) x 4.72 (H); 24 cm (L) x 17 cm (W) x 12cm (H)
Weight	3.96 lbs (1.8 kg)
Filters	Pollen and optional ultra-fine
Device Set-up	On-board keypad and SmartCard [™]
Data Storage Capacity (minimum)	6 month @ 3 usages sessions/day; 7-9 days in-depth data analysis
Data Storage Content	Date/time, duration at pressure, patient event histograms, short-term detailed events, long-term trends, FOSQ tests
Compliance Meter	Breathing detection; sessions > 4 hours
Altitude Compensations	Automatic
Electrical Requirements	100 - 240 AC, 12 VDC
DC Power	Direct-connect cord
Warranty	2 years
Misc.	Integrated heated humidification control, lighted keys on pad, automatic on/off, mask off alert