

RESMED

Share the experience that has helped millions of people like you.

Ultra Mirage II Nasal Mask

Millions of people using positive airway pressure therapy have benefited from ResMed's Ultra Mirage Nasal Mask.

The Ultra Mirage II Nasal Mask combines the best of the Ultra Mirage and other popular ResMed nasal masks.

It offers you a unique combination of comfort, convenience, quietness and effectiveness.



ULTRA MIRAGE™ II
NASAL MASK



Ultra Mirage™ II
Nasal Mask

Patient Information



ResMed Corp Poway, CA, USA +1 858 746 2400 OR 1 800 424 0737 (toll free), ResMed Ltd Bella Vista, NSW, Australia, +61 (2) 8884 1000 or 1 800 658 189 (toll free). Offices in Austria, Finland, France, Germany, Hong Kong, Japan, Malaysia, Netherlands, New Zealand, Singapore, Spain, Sweden, Switzerland, United Kingdom (see website for details). Protected by patents: AU 710733, AU 741003, AU 766623, AU 772832, AU 775051, AU 777033, CA 2261790, DE 29724224, EP 0956069, NZ 513052, NZ 526165, US 6112746, US 6357441, US 6374826, US 6412487, US 6439230, US 6463931, US 6532961, US 6557556, US 6581602, US 6634358, US 6691707, US 6796308. Other patents pending. Protected by design registrations: AU 139764, DE 49911833, DE 40301991, FR 997839, FR 031425, JP 1117921, JP 1197930, US D443355, US D486227, US D493522. Others pending. Ultra Mirage and Mirage are trademarks of ResMed Ltd and Mirage is registered in U.S. Patent and Trademark Office. Specifications may change without notice. ©2005 ResMed Ltd. 1010462/1 05 03

Ultra comfortable

The Ultra Mirage II, with ResMed's Mirage™ dual-wall cushion technology, offers vital comfort.

- The supportive inner layer provides stability.
- The soft outer membrane reduces the likelihood of leaks and pressure sores.
- The cushion helps form a secure seal without the need for tight headgear.

Ultra convenient

- Easy-to-fit, the Ultra Mirage II provides four convenient cushion sizes—Standard, Large, Shallow and Shallow-wide.
- One hand is all it takes to adjust the forehead support and choose one of four positions—numbered and clearly marked for your convenience.

Ultra quiet

One of the quietest masks available, the Ultra Mirage II is designed to help you enjoy restful sleep.

- The built-in air vent directs air away from you and your partner.
- Distinctive design reduces the noise of air released from the vent.

Ultra effective

The Ultra Mirage II gives you extra security and stability.

- The cushion clip keeps the mask cushion firmly attached to the frame.
- Four points of attachment between the mask frame and the headgear enhance stability.


Discover the ResMed difference

ResMed Mirage represents premium, tried and trusted products for positive airway pressure therapy.

The Ultra Mirage has been proving it for over five years—quietly. The Ultra Mirage II continues to redefine what you can expect from the leading nasal mask in the market.



Soft Mirage dual-wall cushion
optimizes seal and minimizes pressure on the bridge of your nose



Adjustable forehead support
offers a choice of four positions for better seal and stability



Flexible, one-piece forehead pads
conform to the shape of your face



Unique air vent design reduces noise by directing the flow of air that is being inhaled and exhaled



360° elbow rotation lets you choose the most convenient tube position

