

## Specifications | REMstar® Auto M Series with C-Flex™

Pressure Range:	4 to 20 cm H <sub>2</sub> O
Ramp Time :	0 to 45 min. (5-min. increments)
Starting Ramp Pressure:	4 to CPAP; patient-adjustable
Dimensions:	7.5" L x 5.0" W x 3.125" H 19 cm L x 12.7 cm W x 7.9 cm H
Weight:	2.2 lbs
Filters:	Pollen and optional ultra-fine
Device Set-up:	LCD/keypad and SmartCard (optional)
Data Storage Capacity:	Display: 7 and 30 day averages SmartCard: 6 months at 3 usage sessions per day, 7 days in-depth analysis
Data Storage Content:	SmartCard: Date/time, hours of use, usage patterns, AHI, leak, snore, FOSQ
Compliance Meter:	Breathing detection
Altitude Compensation:	Automatic
Electrical Requirements:	100 – 240 VAC, 50/60 Hz
Humidification:	Integrated Heated Humidifier and Pass-over Humidifier
DC Power:	Direct connect cord
Warranty:	2 years
Special User Features :	Lighted LEDs, Auto on/off, mask off alert, patient reminders