

Specifications | REMstar[®]

Pressure Range:	4 to 20 cm H ₂ O
Ramp Time :	0 to 45 min. (5-min. increments)
Starting Ramp Pressure:	4 to CPAP; patient-adjustable
Dimensions:	7.5" L x 5.0" W x 3.125" H 19 cm L x 12.7 cm W x 7.9 cm H
Weight:	2.2 lbs
Filters:	Pollen and optional ultra-fine
Device Set-up:	LCD/keypad
Data Storage Capacity (minimum):	Therapy meter
Data Storage Content:	Total therapy hours Total blower hours Usage Sessions > 4 hrs
Compliance Meter:	Breathing detection
Altitude Compensation:	Low, medium or high
Electrical Requirements:	100 – 240 VAC, 50/60 Hz
Humidification:	Integrates with M Series Heated Humidifier and M Series Pass-over Humidifier
DC Power:	Direct connect cord
Warranty:	2 years
Misc:	Fine pressure calibration, lighted LEDs